



T0851PV
SCALETTA A DUE GRADINI

Altezza 42 cm.
Struttura in tubo d'acciaio verniciato con polveri epossidiche. Colore Grigio.
Due gradini in materiale plastico antiscivolo. Colore nero.
Piedini in gomma.



T0851PV
MOUNTING STEP – 2 STEPS

Height 42 cm.
Structure made of epoxy painted steel tube. Grey color.
Two antislip, plastic steps.
Black color.
Rubber feet.



T0851PV
MARCHE-PIED, 2 MARCHES

Hauteur 42 cm.
Structure en tube d'acier verni époxy. Couleur grise.
Deux marches en plastique antidérapante.
Couleur noire.
Pieds en caoutchouc.



T0851PV
ESCABEL, DOS ESCALONES

Altura 42 cm.
Estructura en tubo de acero pintado epoxi. Color gris.
Dos escalones en material plástico antideslizante.
Color negro.
Pies en goma

