



A5-289-120



A5-301-12C



BARATTOLO INOX PORTA COTONE

Interamente stampato in acciaio inox. Coperchio a pomolo.

Mod.	Dim. Esterne
A5-288-100	Ø 100x100 h mm.
A5-289-120	Ø 120x120 h mm.
A5-290-150	Ø 150x150 h mm

BARATTOLO PORTA COTONE
Interamente stampato in acciaio inox. Coperchio a cerniera.

Mod.	Dim. Esterne
A5-300-10C	Ø 100x100 h mm.
A5-301-12C	Ø 120x120 h mm.
A5-302-15C	Ø 150x150 h mm.



S.S. JAR FOR COTTON

Entirely pressed and made of stainless steel. Lid with knob.

Mod.	Overall Dimensions
A5-288-100	Ø 100x100 h mm.
A5-289-120	Ø 120x120 h mm.
A5-290-150	Ø 150x150 h mm

S.S. JAR FOR COTTON
Entirely pressed and made of stainless steel. Hinged lid.

Mod.	Overall Dimensions
A5-300-10C	Ø 100x100 h mm.
A5-301-12C	Ø 120x120 h mm.
A5-302-15C	Ø 150x150 h mm.



BOITE INOX PORTE COTON

Entièrement matricée en acier inoxydable. Couvercle à pommeau

Mod.	Dim. Extérieures
A5-288-100	Ø 100x100 h mm.
A5-289-120	Ø 120x120 h mm.
A5-290-150	Ø 150x150 h mm

BOITE INOX PORTE COTON
Totalemt matricée en acier inoxydable. Couvercle à charnière.

Mod.	Dim. Extérieures
A5-300-10C	Ø 100x100 h mm.
A5-301-12C	Ø 120x120 h mm.
A5-302-15C	Ø 150x150 h mm.



CUBETA INOX PARA ALGODÓN

Totalmente embutido de acero inoxidable. Tapa con bola.

Mod.	Medidas Externas
A5-288-100	Ø 100x100 h mm.
A5-289-120	Ø 120x120 h mm.
A5-290-150	Ø 150x150 h mm

CUBETA INOX PARA ALGODÓN
Totalmente construida de acero inoxidable. Tapa con bisagra.

Mod.	Medidas Externas
A5-300-10C	Ø 100x100 h mm.
A5-301-12C	Ø 120x120 h mm.
A5-302-15C	Ø 150x150 h mm.